

Mobile Hyperbaric System



Innovation & Excellence in Hyperbarics

Mobile Hyperbaric System

PCCI Hyperbaric Systems
300 North Lee Street, #201
Alexandria, VA

Phone: 703-861-8535

Email: info@pccihyperbarics.com

www.pccihyperbarics.com

Applicable Rules

NFPA 99 Chapter 14 "Hyperbaric Facilities"

ASME-PVHO-1 "Safety Code for Pressure Vessels for Human Occupancy"

AIA Guidelines for Mobile Medical Facilities

Chambers


Most Monoplace Chambers Up to 48" Frame Width
or Sized for Customers Needs

Trailer Size

32 to 40 Foot Long

Length Available To Fit Customers Needs

Floor Plans Available To Fit Customers Needs



Easily Transported To Your Site



Easily Sets Up at Your Site



Interior with Two Sechrist 3600 Chambers with PCCI TV System & Gurney



Interior with Two Smaller Sechrist Chambers



Interior with Two PAHI Chambers & Ancillary Items



Interior Showing One Sechrist 3200 Chamber Equipped with PCCI TV & Side Door Wheelchair Lift



Access Options: Doors & Wheelchair Lift



Access Options: Ramp for Gurney Access through Rear Doors Newfoundland, Canada



Interior Options: Countertop with Stainless Steel Sink Hot /Cold Water & Storage Cabinets



Outfitted and Ready for Transport



Initial Operational Location

Hutchinson Hospital, Hutchinson, KS



Treatment in Progress




External Portable LOX Tanks (Used when piped oxygen is not available)



Ready for Treatment Operations San Jose, CA





**For Inquiries
Contact Us Below**

PCCI Hyperbaric Systems
300 North Lee Street, #201
Alexandria, VA
Phone: 703-861-8535
Email: info@pccihyperbarics.com
www.pccihyperbarics.com

Innovation & Excellence in Hyperbarics